Reusable Surgical Instruments

⚠ Important

All products covered by these reprocessing instructions can be found in the appendix.

Before using these products, please read the following information thoroughly!

⚠ Warnings

- Reusable instruments from Ackermann Instrumente GmbH are supplied non-sterile, which is indicated on the product label by the following symbol:

- These products must be cleaned and sterilised before first use in accordance with the instructions in this reprocessing manual.
- The instruments may only be used by trained and, if applicable, instructed physicians. Disinfection, cleaning, and sterilization may only be performed by specially trained medical personal.
- The following instructions only apply to reusable Ackermann products. Please note that any deviation from these instructions, including the use of cleaners / detergents not specifically indicated in these instructions will require an evaluation of device specific efficacy and suitability of the actually performed cycle. Respective evaluation usually requires equipment qualification and device specific performance qualification / validation.
- In case of any serious incident that has occurred in relation to the respective device shall be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established. You can inform us using the form provided on our homepage.
- If product-specific instructions for use are enclosed with the products, these must also be followed.

Intended Use

Reusable surgical instruments manufactured by Ackermann Instrument GmbH are intended to perform general surgical procedures. It is the responsibility of the attending physician to determine the use of each product for a specific operation. The intended use of the respective product group can be found in the appendix.

Contraindication(s)

The instruments are contraindicated for use on the central circulatory system and central nervous system.

Not for use in patients allergic to NI, CR, brass or aluminum.

The products are contraindicated for all applications other than the intended use described.

If there are any product-specific contraindications, these are listed in the instructions for use of the respective product.

Limitations on Reprocessing

The reusable surgical instruments are designed for repeated use. The service life of the instruments results from reprocessing and handling. Therefore, a limit for the number of reprocessing cycles cannot be specified. Before each use, the product must be checked for function and damage. Damaged products must be disposed of.

⚠ Warning

If the instruments are used on patients infected with prion-based diseases (e.g. CJD) or HIV, the instrument must not be reused. Even after reprocessing and sterilization, the risk of cross-contamination cannot be eliminated and the instrument must be destroyed.

Inspection before Use

⚠ Warning

Before operating the instruments, a visual inspection and a functional test must be carried out in accordance with section "Functional control and maintenance of the products in the reprocessing cycle according to DIN 96298-4". Non-functional instruments must be discarded.

Reprocessing Instructions

⚠ Warning

Always ensure that the products are used and reprocessed by qualified personnel with appropriate experience in hospital hygiene and sterilization techniques. To ensure safe and effective reprocessing of the products, the following instructions have been validated by the manufacturer for effectiveness and compatibility with the instruments. It is the responsibility of the end user to ensure that cleaning and sterilization is indeed performed with appropriate equipment, materials and personnel to achieve the desired result. Any deviation from these instructions should be evaluated for effectiveness and possible adverse consequences.

⚠ Warning

All reusable Ackermann products must be subjected to reprocessing as described in the following sections before being used for the first time. Follow the instructions and warnings of the manufacturers of the decontamination, disinfection and cleaning agents to be used.

Preparation before Use

⚠ Warning

Ackermann reusable products are supplied non-sterile and must be cleaned and sterilized before first use and before each subsequent use. The packaging cannot withstand the high temperatures of autoclaving and should be disposed of before sterilization.

Cleaning

⚠ Warning

Improper cleaning, rinsing, and drying can result in potentially hazardous residues or inadequate sterilization.

Pre-Cleaning

Instruments, or fully disassembled instruments if applicable, must be brushed under cold water until all visible contamination is removed. If the instrument has a flushing port, the lumen of the shaft must be flushed for at least 10 seconds through the flushing port with a water jet gun (static pressure greater than 4 bar) after manual brushing. The instruments are then treated in a 0.5% alkaline enzymatic (e.g. Neodisher® Mediclean forte) ultrasonic bath for 10 minutes. Finally, rinse again for 30 seconds with a water jet gun.

Automated Cleaning¹

The automated cleaning was validated with Neodisher® Mediclean forte (alkaline-enzymatic: pH-value 10,4 – 10,8) for cleaning and Neodisher® Z (acidic: pH-value 3,0 – 2,6) for neutralization.

| Step | Parameter | |
|-----------|----------------------|--|
| Pre-rinse | Rinsing temperature | Cold tap water |
| | Exposure time | 180 secs. / 3 mins. |
| Cleaning | Cleaning temperature | 55°C |
| | Exposure time | 600 secs. / 10 mins. |
| | Detergent | Alkaline-enzymatic (e.g. Neodisher® Mediclean forte) |
| | Concentration | 0.50% |

| | | |
|----------------------|--------------------------|----------------------------|
| | Rinsing temperature | Cold DI water |
| | Exposure time | 60 secs. / 1 min. |
| | Neutralizing agent | acidic (e.g. Neodisher® Z) |
| | Concentration | 0.10% |
| Rinse | Rinsing temperature | Cold DI water |
| | Exposure time | 120 secs. / 2 mins. |
| Thermal disinfection | Disinfection temperature | 90 °C (A0 3000) |
| | Exposure time | 300 secs. / 5 mins. |
| Drying | Drying | 60 mins. / 1 hour |

Steam Sterilization²

| | |
|-------------|--------------------------------|
| Procedure | Triple fractionated pre-vacuum |
| Temperature | min. 132°C |
| Duration | min. 3 mins. |
| Drying time | min. 20 mins. |

The instructions provided above have been validated as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

¹ Validated with washer-disinfector Miele PG 8535

² Validated with steam autoclave Lautenschläger ZentraCert

Function Control and Maintenance of the Products in the Reprocessing Cycle acc. to DIN 96298-4

⚠ Warning

If the product is damaged or suspected of being damaged, do not attempt to repair the product yourself. Do not reuse damaged instruments.

⚠ Warning

If a product-specific inspection or maintenance is required, the relevant information can be found in the instruction of use for the respective product.

Maintenance

Maintenance is always performed prior to functional testing. Oils can be applied by squeeze, spray, or pipe. Oil the joints and sliding surfaces in a directed way, then move the joint and sliding surfaces several

times and wipe off excess oil with a lint-free cloth.

Visual and Fuctional Inspection

Products marked as nonfunctional must be discarded. The following tests and maintenance must be performed by qualified personnel.

Visual Inspection

Visual inspection for cleanliness

- The instruments should be free of macroscopic contaminants.

Visual control for missing parts

- All parts such as screws, pins, springs, valves, hooks, etc. must be present. If this is not the case, the whereabouts of the missing parts must be clarified.

Visual control for breakage

- If parts (such as tips, teeth, springs, working ends, TC inserts) are missing due to material breakage, the whereabouts of the missing parts must be clarified.

Check for deformation

- Instruments with deformed working ends, cutting edges, tips, etc. need to be sorted out.

Visual inspection for corrosion

- Instruments with visible corrosion in the form of pitting, fretting, stress or crevice corrosion must be sorted out.

Visual inspection for discoloration

- Instruments with a discoloration must be sorted out.

Visual inspection for cracks

- Instruments with a crack must be sorted out.

Visual inspection of coatings

- Coatings must be free of imperfections and notches. If damage to a coating is detected, this instrument must be sorted out.

Visual inspection of markings

- The marking must be legible. If the instrument is marked with a machine-readable marking, this must be readable with an appropriate scanner. If the marking is no longer legible, the instrument must be sorted out.

Fuctional Check

Function check for clamps

- The clamps must have a uniform, backlash-free movement.
- The lock must hold securely in any position.
- The working ends must be congruent with each other.
- The teeth, if any, must mesh without hooking.

Function check for scissors

- When closing the scissors, there must be minimal and continuous friction between the scissor blades, but no noticeable snagging. A cutting test can be carried out to eliminate any doubts.

Function check for needleholders

- The needleholder must have a smooth, backlash-free movement.
- The catch must hold securely in any position.
- The working ends must be congruent.
- The working ends must close completely; if a gap is detected, the tooth profile must be checked for wear and replaced if necessary.

Function check for tweezers

- The teeth must engage without hooking. This is checked by pressing the tweezer legs together. The tips of the forceps (grasping surfaces) must be congruent with each other and must not slide sideways from each other when the forceps legs are pressed together.

Function check for retractors

- The retractor must hold securely in any position during opening. This is checked by opening and simultaneously applying pressure to the arms. Valves must be easy to insert and still fit securely.

Functional check for chisels and osteotomes with handles made of plastics or aluminum

- The components of the instrument must be firmly connected to each other without wobbling. The working ends must be checked for the presence of nicks. Impact caps with severe deformation or splintering must be discarded.

Function check for bone splitters and rongeurs

- The rongeurs with spring must open independently. The cutting edges should not have any major nicks, otherwise the function according to the intended purpose is no longer guaranteed. Rongeurs with severely damaged cutting edges must be discarded.

Labeling and Storage

Labeling

⚠ Warning

Please pay attention to the symbols on the product label and in this instruction for use.

Care must be taken to ensure that the instruments are fully traceable at all times.

Storage Conditions

Storage conditions can be found on the product label.

⚠ Warning

The storage area must be dust-free, low in microbiological contamination, dark and free from temperature fluctuations.

Disposal

Environmentally sound disposal enables valuable raw materials to be recycled. Dispose of the device in an environmentally friendly manner in accordance with the valid hospital guidelines.

Repairs

Despite being used as intended, medical devices are subject to varying degrees of wear depending on the intensity of use.

- Please note that repairs may only be carried out by the manufacturer or authorised personnel.
- The medical devices must be cleaned, disinfected and sterilised before being sent in for repair. Soiled or contaminated medical devices must not be sent.

Additional Information

Do not exceed maximum loading capacity of the sterilizer when processing multiple instruments in one sterilization cycle.

Packaging

Packaging suitable for steam sterilization must comply with the requirements according to ISO 11607 / ANSI / AAMI ST79 / AAMI TIR12, for example, disposable sterilization packs (single or double packs) temperature resistant up to at least 137°C (279°F) and sufficient steam permeability, which provide sufficient protection against mechanical damage, or sterilization containers which need to be maintained according to the manufacturer's instructions.

Warranty

A guarantee of two years from handover to the end customer is given for material and manufacturing defects. Transport costs and shipping risk are not assumed. In addition, the warranty stated in the General Terms and Conditions applies. This warranty does not apply to instruments that have been repaired, modified or improperly handled by unauthorised persons. In this case, no liability is assumed for the operational safety of the instrument either. In addition, the warranty becomes null and void.

⚠ Warning

Risk of infection due to non-sterile instruments.

- Prepare the instrument before returning it to the manufacturer.
- Return the reconditioned instrument in its original packaging to the manufacturer.

Symbols used on labeling (acc. ISO 15223-1)

| | |
|---|--|
|  | Legal manufacturer |
| LOT | Batch code |
| REF | Product number |
|  | Non-sterile |
|  | Do not use if package is damaged |
|  | Keep dry |
| MD | Medical device |
| UDI | Unique device identifier |
|  | Consult instruction for use |
|  | Caution |
|  | Keep away from sunlight |
| R^{ONLY} | Only by prescription from medical doctor |
| QTY | Quantity |
| CE 0483 | mdc medical device certification GmbH Kriegerstraße 6 70191 Stuttgart, Germany |

Contact

Ackermann Instrumente GmbH 🏢
 Eisenbahnstrasse 65-67
 78604 Rietheim-Weilheim
 Germany

Phone: +49 (0)7461 966 17 - 0
 Fax: +49 (0)7461 966 17 - 70
 E-mail: sales@ackermannsurgical.com
 Web: www.ackermannsurgical.com

End of document

Further Information

ROLAN is a trademark of Ackermann Instrumente GmbH.

| |
|----------|
| Appendix |
|----------|

| Product Group | Intended Use | Basic UDI |
|--|---|-----------------|
| Cannula, reusable, injection, non-sterile, laparoscopy | The instrument is intended to be manually-operated to perform various needle functions (e.g., injection, aspiration, irrigation, puncture). | 42503022034A0N4 |
| Guiding instrument, reusable, non-sterile, arthroscopy | The instrument is intended to simultaneously place, angle, and guide a rotating surgical drill on hard tissue (e.g., bone). | 42503022045A0NG |
| Endoscopic forceps, flexible, reusable, biopsy, non-sterile, gynecology | The instrument is intended to obtain tissue specimens for examination to determine the presence or extent of a disease. | 42503022046A0NM |
| Endoscopic forceps, rigid, reusable biopsy, non-sterile, arthroscopy | The instrument is intended to obtain tissue specimens for examination to determine the presence or extent of a disease. | 42503022047A0NS |
| Endoscopic forceps, rigid, reusable biopsy, non-sterile, bronchoscopy | The instrument is intended to obtain tissue specimens for examination to determine the presence or extent of a disease. | 42503022048A0NX |
| Endoscopic forceps, rigid, reusable biopsy, non-sterile, laparoscopy | The instrument is intended to obtain tissue specimens for examination to determine the presence or extent of a disease. | 42503022049A0P4 |
| Endoscopic forceps, rigid, reusable biopsy, non-sterile, urology | The instrument is intended to obtain tissue specimens for examination to determine the presence or extent of a disease. | 42503022050A0MW |
| Endoscopic forceps, rigid, reusable cutting/punching, non-sterile, ENT | The instrument is intended to punch/cut/excise soft-tissue. | 42503022051A0N3 |
| Endoscopic forceps, rigid, reusable cutting/punching, non-sterile, arthroscopy | The instrument is intended to punch/cut/excise soft-tissue. | 42503022052A0N8 |
| Dilator, reusable, non-sterile, arthroscopy | The instrument is intended to dilate a body lumen or access channel. | 42503022053A0ND |
| Drilling instrument, reusable, non-sterile, arthroscopy | The instrument is intended to cut and remove circular bone sections by rotary motions. | 42503022054A0NJ |
| Endoscopic retractor, reusable, non-sterile, laparoscopy | The instrument is intended to lift urogenital structures. | 42503022055A0NP |
| Endoscopic forceps, flexible, reusable, grasping, non-sterile, gynecology | The instrument is intended to grasp and manipulate tissue or additionally foreign bodies. | 42503022056A0NU |
| Endoscopic forceps, flexible, reusable, grasping, non-sterile, urology | The instrument is intended to grasp and manipulate tissue or additionally foreign bodies. | 42503022057A0NZ |
| Endoscopic forceps, flexible, reusable, grasping, non-sterile, arthroscopy | The instrument is intended to grasp and manipulate tissue or additionally foreign bodies. | 42503022058A0P6 |
| Endoscopic forceps, rigid, reusable grasping, non-sterile, bronchoscopy | The instrument is intended to grasp and manipulate tissue or additionally foreign bodies. | 42503022059A0PB |
| Endoscopic forceps, rigid, reusable grasping, non-sterile, ENT | The instrument is intended to grasp and manipulate tissue or additionally foreign bodies. | 42503022060A0N5 |
| Endoscopic forceps, rigid, reusable grasping, non-sterile, laparoscopy | The instrument is intended to grasp and manipulate tissue or additionally foreign bodies. | 42503022061A0NA |
| Endoscopic forceps, rigid, reusable grasping, non-sterile, urology | The instrument is intended to grasp and manipulate tissue or additionally foreign bodies. | 42503022062A0NF |
| Impact instrument, reusable, non-sterile, arthroscopy | The instrument is intended to position an orthopaedic implant within the body by transmission of an impact force that drives the implant into its final position. | 42503022063A0NL |
| Knife, reusable, non-sterile, arthroscopy | The Instrument is intended to manually cut or dissect tissue. | 42503022066A0P3 |
| Knife, reusable, non-sterile, ENT | The Instrument is intended to manually cut or dissect tissue. | 42503022067A0P8 |
| Knife, reusable, urethrotome system, non-sterile, urology | The Instrument is intended to manually cut or dissect tissue. | 42503022068A0PD |
| Endoscopic forceps, rigid, reusable optical forceps, non-sterile, bronchoscopy | The instrument is intended to grasp and manipulate tissue or additionally foreign bodies. | 42503022069A0PJ |
| Probe, reusable, non-sterile, arthroscopy | The instrument is intended to explore, measure, and/or perform therapeutic procedures. | 42503022070A0NC |
| Ring applicator, reusable, non-sterile, laparoscopy | The instrument is intended to apply small elastic rings to block the fallopian tubes. | 42503022071A0NH |
| Endoscopic scissors, flexible, reusable, non-sterile, gynecology | The instrument is intended to cut tissue or sutures. | 42503022072A0NN |
| Endoscopic scissors, flexible, reusable, non-sterile, urology | The instrument is intended to cut tissue or sutures. | 42503022073A0NT |
| Endoscopic scissors, rigid, reusable, non-sterile, urology | The instrument is intended to cut tissue or sutures. | 42503022074A0NY |
| Endoscopic scissors, flexible, reusable, non-sterile, Bronchoscopy | The instrument is intended to cut tissue or sutures. | 42503022075A0P5 |
| Endoscopic scissors, rigid, reusable non-sterile, laparoscopy | The instrument is intended to cut tissue or sutures. | 42503022076A0PA |
| Hammer, reusable, surgical, non-sterile, arthroscopy | The Instrument is intended to strike a surface or another device. | 42503022077A0PF |
| Surgical procedure kit, reusable, non-sterile, arthroscopy | The instrument kit is intended to perform a specific surgery. | 42503022078A0PL |
| Screwdriver, reusable, non-sterile, arthroscopy | The instrument is intended to fit into a screw head to introduce/remove a screw. | 42503022079A0PR |
| Endoscopic forceps, rigid, reusable suture passer, non-sterile, arthroscopy | The instrument is intended to pass sutures or a needle/suture assembly through tissues/prosthetic materials to facilitate suturing. | 42503022080A0NK |
| Endoscopic forceps, rigid, reusable suture passer, non-sterile, laparoscopy | The instrument is intended to pass sutures or a needle/suture assembly through tissues/prosthetic materials to facilitate suturing. | 42503022081A0NQ |
| Knot pusher, reusable, non-sterile, arthroscopy | The instrument is intended to advance extracorporeal suture knots to the point of appropriate suture tension. | 42503022082A0NV |
| Knot pusher, reusable, non-sterile, laparoscopy | The instrument is intended to advance extracorporeal suture knots to the point of appropriate suture tension. | 42503022083A0P2 |
| Myoma screw, reusable, non-sterile, gynecology | The instrument is intended to engage, stabilize, and manipulate a myoma (i.e., a uterine fibroid) and other gynaecological tissues during their surgical removal. | 42503022084A0P7 |
| Awl, reusable, non-sterile, arthroscopy | The instrument is intended to bore holes in bone. | 42503022085A0PC |
| Curette, reusable, non-sterile, arthroscopy | The Instrument is intended for scraping out or off bones. | 42503022086A0PH |
| Curette, reusable, non-sterile, ENT | The Instrument is intended for scraping out or off bones. | 42503022087A0PN |
| Tendon stripper, reusable, non-sterile, arthroscopy | The Instrument is intended to excise a length of ligament, tendon or fascia for use as living graft. | 42503022088A0PT |
| Rasp, reusable, bone, non-sterile, arthroscopy | The Instrument is intended to be used to scrape and smooth the surface of bone. | 42503022090A0NS |
| Endoscopic forceps, rigid, reusable clip applier, non-sterile, ligation system | The instrument is intended to apply small atraumatic clips for the ligation of blood vessels, or similar tubular anatomy. | 42503022091A0NX |
| Forceps, open-surgery, reusable, clip applier, non-sterile, ligation system | The instrument is intended to apply small atraumatic clips for the ligation of blood vessels, or similar tubular anatomy. | 42503022092A0P4 |

| | | |
|--|--|-----------------|
| Endoscopic forceps, rigid, reusable clip remover, non-sterile, ligation system | The instrument is intended to grasp, open, and remove a tissue clip. | 42503022093A0P9 |
| Forceps, open surgery, reusable, clip remover, non-sterile, ligation system | The instrument is intended to grasp, open, and remove a tissue clip. | 42503022094A0PE |
| Probe, reusable, non-sterile, laparoscopy | The instrument is intended to explore, measure, and/or perform therapeutic procedures. | 42503022095A0PK |
| Retractor, reusable, hand-held, non-sterile, laparoscopy | The instrument is intended to separate/draw aside the margins of a wound/incision to allow access to tissues/organs. | 42503022096A0PQ |
| Retractor, reusable, system component, non-sterile, retraction and positioning | The Instrument is intended to separate the margins of a wound/incision. | 42503022097A0PV |
| Elevator/lever, reusable, non-sterile, ENT | The instrument is intended to lift, position or pry anatomical structures or surgical material. | 42503022098A0Q2 |
| Syringe, reusable, bladder, non-sterile, urology | The product is intended for the injection or withdrawal of liquids or gases into or out of the body. | 42503022099A0Q7 |
| Syringe adapter, reusable, non-sterile, urology | The adapter is intended to enable the interconnection of compatible components. | 42503022099A1Q9 |
| Bladder evacuator, reusable, non-sterile, urology | The Product is intended to suction tissue and body fluids. | 42503022102A0MG |
| Endoscopic forceps, rigid, reusable grasping, non-sterile, VATS | The instrument is intended to grasp and manipulate tissue or additionally foreign bodies. | 42503022112A0MP |
| Endoscopic scissors, rigid, reusable non-sterile, VATS | The instrument is intended to cut tissue or sutures. | 42503022113A0MU |
| Dilator, reusable, non-sterile, urology | The instrument is intended to dilate a body lumen or access channel. | 42503022116A0NB |
| Knot pusher, reusable, non-sterile, VATS | The instrument is intended to advance extracorporeal suture knots to the point of appropriate suture tension. | 42503022118A0NM |
| Endoscopic forceps, rigid, reusable biopsy, non-sterile, VATS | The instrument is intended to obtain tissue specimens for examination to determine the presence or extent of a disease. | 42503022119A0NS |
| Endoscopic forceps, rigid, reusable needle holder, non-sterile, laparoscopy | The Instrument is intended to grasp a suture needle when pushing/pulling the needle and attached suture through tissue during suturing. | 42503022120A0ML |
| Endoscopic forceps, rigid, reusable needle holder, non-sterile, VATS | The Instrument is intended to grasp a suture needle when pushing/pulling the needle and attached suture through tissue during suturing. | 42503022121A0MR |
| Suction/irrigation cannula, rigid, reusable non-sterile, VATS | The instrument is intended for irrigation and/or aspiration of body cavities/wounds. | 42503022122A0MW |
| Endoscopic forceps, rigid, reusable suture cutter, non-sterile, VATS | The instrument is intended for cutting suture. | 42503022123A0N3 |
| Speculum, reusable, nasal, non-sterile, ENT | The instrument is intended to create a channel for examination, suction, irrigation, or the insertion of another surgical device in a body opening. | 42503022130A0MT |
| Knife, reusable, non-sterile, urology | The Instrument is intended to manually cut or dissect tissue. | 42503022131A0MY |
| Endoscopic forceps, flexible, reusable, biopsy, non-sterile, urology | The instrument is intended to obtain tissue specimens for examination to determine the presence or extent of a disease. | 42503022133A0NA |
| Gall bladder extractor, reusable, non-sterile, laparoscopy | The instrument is intended to remove the gallbladder. | 42503022134A0NF |
| Clamp, reusable, penis clamp, non-sterile, urology | The instrument is intended to grasp the penis for manipulation during surgery. | 42503022135A0NL |
| Suture cannula, reusable, non-sterile, laparoscopy | The instrument is intended to pass sutures or a needle/suture assembly through tissues/prosthetic materials to facilitate suturing. | 42503022136A0NR |
| Awl, reusable, non-sterile, general surgery | The instrument is intended to bore holes in bone. | 42503022137A0NW |
| Drill, reusable, orthopaedics, non-sterile, general surgery | The instrument is intended to cut and remove circular bone sections by rotary motions. | 42503022138A0P3 |
| Drilling instrument, reusable, non-sterile, general surgery | The instrument is intended to cut and remove circular bone sections by rotary motions. | 42503022139A0P8 |
| Hammer, reusable, surgical, non-sterile, general surgery | The Instrument is intended to strike a surface or another device. | 42503022140A0N2 |
| Dermatome, reusable, non-sterile, general surgery | The Instrument is intended to manually cut or dissect tissue. | 42503022141A0N7 |
| Dilator, reusable, bile, non-sterile, general surgery | The instrument is intended to dilate a body lumen or access channel. | 42503022142A0NC |
| Dilator, reusable, vascular, non-sterile, general surgery | The instrument is intended to dilate a body lumen or access channel. | 42503022143A0NH |
| Dissector / ENUCLEATOR, reusable, non-sterile, general surgery | The instrument is intended to separate a soft tissue or body structure from another. | 42503022144A0NN |
| Elevator/lever, reusable, ENT, non-sterile, general surgery | The instrument is intended to lift, position or pry anatomical structures or surgical material. | 42503022145A0NT |
| Elevator/lever, reusable, non-sterile, general surgery | The instrument is intended to lift, position or pry anatomical structures or surgical material. | 42503022146A0NY |
| Guide, reusable, suture, non-sterile, general surgery | The Instrument is intended for passing a length of a flexible material/device, though a surgical site. | 42503022147A0P5 |
| File, reusable, orthopaedic, non-sterile, general surgery | The Instrument is intended to scrape and smooth the surface of bone. | 42503022148A0PA |
| Vascular approximator, reusable, non-sterile, general surgery | The instrument is intended to grasp the torn or resected portions of a vessel to hold the portions in apposition for the purpose of reattaching one portion to the other in order to repair the structure. | 42503022149A0PF |
| Implantation instrument, reusable, non-sterile, general surgery | The instrument is intended to aid navigation/placement of an implant/prosthesis or other instruments during implantation. | 42503022150A0N9 |
| Cannula, reusable, injection, non-sterile, general surgery | The instrument is intended to be inserted beneath the skin of patient while connected to a noninvasive device (e.g., syringe, secondary medication set) to administer and/or withdraw (aspirate) fluids/drugs. | 42503022151A0NE |
| Intubation tube, reusable, non-sterile, general surgery | The instrument is intended to be inserted into an artificial opening made in the trachea (windpipe) during tracheotomy to provide a patent airway. | 42503022152A0NK |
| Clamp, reusable, arterial, non-sterile, general surgery | The instrument is intended to grasp, manipulate, dissect, and/or clamp soft-tissues. | 42503022153A0NQ |
| Clamp, reusable, clip applier, non-sterile, general surgery | The instrument is intended to apply small atraumatic clips for the ligation of blood vessels, or similar tubular anatomy. | 42503022154A0NV |
| Clamp, reusable, vascular, non-sterile, general surgery | The instrument is intended to grasp, manipulate, dissect, and/or clamp soft-tissues. | 42503022155A0P2 |
| Clamp, reusable, grasping, non-sterile, general surgery | The instrument is intended to grasp and manipulate tissue or additionally foreign bodies. | 42503022156A0P7 |
| Clamp, reusable, penis clamp, non-sterile, general surgery | The instrument is intended to grasp the penis for manipulation during surgery. | 42503022157A0PC |
| Clamp, reusable, sponge/bandage/cloth, non-sterile, general surgery | The instrument is intended to be used for non-dedicated grasping of devices, sponges and/or dressings during a procedure. | 42503022158A0PH |
| Clamp, reusable, bone holding, non-sterile, general surgery | The instrument is intended to grasp and hold a bone. | 42503022159A0PN |
| Curette, reusable, non-sterile, general surgery | The Instrument is intended for scraping out or off bones. | 42503022160A0NG |
| Marker, reusable, corneal, non-sterile, general surgery | The instrument is intended to imprint, indent, and/or incise corneal tissue prior to an ophthalmic surgical procedure. | 42503022161A0NM |

| | | |
|---|--|-----------------|
| Chisel, reusable, non-sterile, general surgery | The instrument is intended to remove or cut through bone. | 42503022162A0NS |
| Knife, reusable, non-sterile, general surgery | The Instrument is intended to manually cut or dissect tissue. | 42503022163A0NX |
| Clamp, reusable, umbilical cord, non-sterile, general surgery | The instrument is intended to temporarily compress the umbilical cord immediately after birth. | 42503022165A0P9 |
| Needle, reusable, general-purpose, non-sterile, general surgery | The instrument is intended to be manually-operated to perform various needle functions (e.g., injection, aspiration, irrigation, puncture). | 42503022166A0PE |
| Guide needle, reusable, non-sterile, general surgery | The instrument is intended to drain fluids. | 42503022166A0PE |
| Clamp, reusable, needle holder, non-sterile, general surgery | The Instrument is intended to grasp a suture needle when pushing/pulling the needle and attached suture through tissue during suturing. | 42503022168A0PQ |
| Needle, reusable, suture, non-sterile, general surgery | The instrument is intended to insert and pull suture material through tissue. | 42503022169A0PV |
| Osteotome, reusable, non-sterile, general surgery | The Instrument is intended to cut and/or shape bone by percussion during an orthopaedic procedure. | 42503022170A0NP |
| Tweezers, reusable, grasping, non-sterile, general surgery | The instrument is intended to grasp and manipulate tissue or additionally foreign bodies. | 42503022171A0NU |
| Tweezers, reusable, grasping, ophthalmology, non-sterile, general surgery | The instrument is intended to grasp and manipulate tissue or additionally foreign bodies. | 42503022172A0NZ |
| Rasp, reusable, bone, non-sterile, general surgery | The instrument is intended to scrape and smooth the surface of bone. | 42503022173A0P6 |
| Rasp, reusable, nose, non-sterile, general surgery | The instrument is intended to remove and profile the cartilage/bone of the nose. | 42503022174A0PB |
| Retractor, reusable, self-retaining, non-sterile, general surgery | The instrument is intended to separate/draw aside the margins of a wound/incision to allow access to tissues/organs. | 42503022175A0PG |
| Retractor, reusable, hand-held, non-sterile, general surgery | The instrument is intended to separate/draw aside the margins of a wound/incision to allow access to tissues/organs. | 42503022176A0PM |
| Retractor, reusable, ophthalmology, hand-held, non-sterile, general surgery | The instrument is intended to separate/draw aside the margins of a wound/incision to allow access to tissues/organs. | 42503022178A0PX |
| Retractor, reusable, system component, self-retaining, non-sterile, general surgery | The Instrument is intended to separate the margins of a wound/incision. | 42503022179A0Q4 |
| Retractor, reusable, ophthalmology, self-retaining, non-sterile, general surgery | The instrument is intended to separate/draw aside the margins of a wound/incision to allow access to tissues/organs. | 42503022180A0NW |
| Rongeur, reusable, non-sterile, general surgery | The instrument is intended for cutting tough tissues (intervertebral disk) or bone, through a cutting/biting action. | 42503022181A0P3 |
| Saw, reusable surgical, manual, non-sterile, general surgery | The instrument is intended for cutting bones or other body tissue by means of a sawing action performed by the user. | 42503022182A0P8 |
| Saw, reusable, nose, non-sterile, general surgery | The instrument is intended for cutting bone and/or cartilage in the nose. | 42503022183A0PD |
| Suction/irrigation cannula, rigid, reusable non-sterile, general surgery | The instrument is intended for irrigation and/or aspiration of body cavities/wounds. | 42503022184A0PJ |
| Scissors, open-surgery, reusable, non-sterile, general surgery | The instrument is intended to cut tissue or sutures. | 42503022185A0PP |
| Snare, reusable, nasal, non-sterile, general surgery | The instrument is intended to be inserted into the naris for the removal of tissue, typically polyps, tumours, and other abnormal tissue from the nasal cavity during ear/nose/throat (ENT) surgery. | 42503022186A0PU |
| Snare, reusable, ear, non-sterile, general surgery | The instrument is intended to be inserted through the outer ear canal or through a surgical incision in the ear to remove tissue, typically tumorous or damaged tissue, from the ear during ear/nose/throat (ENT) surgery. | 42503022187A0PZ |
| Snare, reusable, tonsil, non-sterile, general surgery | The Instrument is intended to remove the tonsils during ear/nose/throat (ENT) surgery. | 42503022188A0Q6 |
| Screwdriver, reusable, non-sterile, general surgery | The instrument is intended to fit into a screw head to introduce/remove a screw. | 42503022189A0QB |
| Tendon stripper, reusable, non-sterile, general surgery | The Instrument is intended to excise a length of ligament, tendon or fascia for use as living graft. | 42503022190A0P5 |
| Probe, reusable, non-sterile, general surgery | The instrument is intended to explore, measure, and/or perform therapeutic procedures. | 42503022191A0PA |
| Spatula/needle, reusable, ophthalmology, non-sterile, general surgery | The instrument is intended to remove a foreign body/object embedded in or adhering to the surface of the eye globe. | 42503022192A0PF |
| Spatula, open-surgery, reusable, non-sterile, general surgery | The instrument is intended to manipulate tissue during various types of surgical interventions. | 42503022193A0PL |
| Forceps, open-surgery, reusable, skin punching-, general surgery | The instrument is intended to punch/cut/excise soft-tissue. | 42503022194A0PR |
| Syndesmotome, reusable, non-sterile, general surgery | The instrument is intended for cutting desmodontal fibres and for detaching the alveolodental ligaments. | 42503022195A0PW |
| Marker, reusable, ophthalmologic, non-sterile, general surgery | The instrument is intended to measure the diameter, length, angles, and thicknesses of eye. This instrument is not calibrated. | 42503022196A0Q3 |
| Trephine, reusable, corneal, non-sterile, general surgery | The instrument is intended to cut and remove circular pieces of corneal tissue. | 42503022197A0Q8 |
| Trocar, reusable, general-purpose, non-sterile, general surgery | The instrument is intended to create a body cavity by puncture to facilitate the introduction of an access cannula. | 42503022198A0QD |
| Trocar, reusable, tracheal, non-sterile, general surgery | The instrument is intended to create a body cavity by puncture to facilitate the introduction of an access cannula. | 42503022199A0QJ |
| Catheter, reusable, ear, non-sterile, general surgery | The instrument is intended for middle ear drainage and/or temporary equalization of the pressure on both sides of the eardrum via insufflation. | 42503022200A0MH |
| Forceps, open surgery, reusable, biopsy, non-sterile, general surgery | The instrument is intended to obtain tissue specimens for examination to determine the presence or extent of a disease. | 42503022201A0MN |
| Forceps, open surgery, reusable, punch, ENT, non-sterile, general surgery | The instrument is intended to punch/cut/excise soft-tissue. | 42503022202A0MT |
| Forceps, open surgery, reusable, grasping, non-sterile, general surgery | The instrument is intended to grasp and manipulate tissue or additionally foreign bodies. | 42503022203A0MY |
| Forceps, open surgery, reusable, bone holding, non-sterile, general surgery | The instrument is intended to grasp and hold a bone. | 42503022204A0N5 |
| Forceps, open surgery, reusable, bone punch, non-sterile, general surgery | The instrument is intended to punch/cut/excise posterior rib ends and bone tissue. | 42503022205A0NA |
| Forceps, reusable, tooth extraction, non-sterile, general surgery | The instrument is intended for the extraction of teeth. | 42503022206A0NF |
| Forceps, open surgery, reusable, wire forceps, non-sterile, general surgery | The instrument is intended to cut, grip, tighten and/or twist wires that are being applied to the patient. | 42503022207A0NL |
| Intubation laryngoscope rigid, reusable non-sterile, general surgery | The instrument is intended to manipulate the tongue, preventing it from obstructing the oropharynx and enabling a clear view of the trachea. | 42503022216A0NN |
| Speculum, reusable, nasal, non-sterile, general surgery | The instrument is intended to create a channel for examination, suction, irrigation, or the insertion of another surgical device in a body opening. | 42503022218A0NY |
| Speculum, vaginal, reusable, general surgery | The instrument is intended to create a channel for examination, suction, irrigation, or the insertion of another surgical device in a body opening. | 42503022219A0P5 |
| Dilator, reusable, urethral, non-sterile, general surgery | The instrument is intended to dilate a body lumen or access channel. | 42503022220A0MX |

| | | |
|--|--|-----------------|
| Speculum, reusable, ear, non-sterile, general surgery | The instrument is intended to create a channel for examination, suction, irrigation, or the insertion of another surgical device in a body opening. | 42503022221A0N4 |
| Mouth gag, reusable, non-sterile, general surgery | The instrument is intended to be placed between the teeth of the upper and lower jaws of a patient to maintain an open oral cavity. | 42503022222A0N9 |
| Catheter, reusable, ureteral, non-sterile, general surgery | The device is intended for ureter dilation, bypass of partial obstructions, urinary tract irrigation/drainage during surgery, stone removal, retrograde urography, guidewire placement, and/or brush biopsy. | 42503022223A0NE |
| Dilator, uterine, reusable, general surgery | The instrument is intended to dilate a body lumen or access channel. | 42503022224A0NK |
| Dental hand instrument, reusable, non-sterile, general surgery | The dental hand instruments are intended in dentistry to perform various diagnostic, preventive and therapeutic procedures. | 42503022225A0NQ |
| Dilator, reusable, lacrimal, non-sterile, general surgery | The instrument is intended to dilate a body lumen or access channel. | 42503022228A0P7 |
| Dilator, reusable, tracheal surgery, non-sterile, general surgery | The instrument is intended to dilate a body lumen or access channel. | 42503022229A0PC |
| Speculum, reusable, rectal, non-sterile, general surgery | The instrument is intended to create a channel for examination, suction, irrigation, or the insertion of another surgical device in a body opening. | 42503022234A0NS |
| Syringe, reusable, general-purpose, non-sterile, general surgery | The product is intended for the injection or withdrawal of liquids or gases into or out of the body. | 42503022235A0NX |
| Endoscopic forceps, rigid, reusable cutting, ENT, non-sterile, general surgery | The instrument is intended to be used to punch/cut/excise soft-tissue. | 42503022237A0P9 |
| Surgical procedure kit, reusable, non-sterile, general surgery | The instrument kit is intended to perform a specific surgery. | 42503022239A0PK |
| Optical tube, reusable, rectal, non-sterile, general surgery | The cannula is intended to provide access from a small body opening (surgical/non-surgical) to the surgical site to allow access for endoscopes, instruments and/or fluids. | 42503022240A0ND |
| Spatula, non-surgical, reusable, non-sterile, general surgery | The instrument is intended to manipulate tissue. The instrument is intended to be non-surgical. | 42503022253A0P3 |

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CE mdc medical device certification GmbH
ROLAN ist eine Marke der Ackermann Instrumente GmbH.